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Chapter: Workflow (All)

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Chapter: System Level Checks (Visits 1,3)

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Engineer Name: Takuya Nishio

Device Status: Pass

Event Type: Predictive / Preventative Maintenance

Service Activity: Visit 4

EQUIPMENT DETAILS

Product Name: Ingenia 3.0T

Product Number: 781377

Serial Number: 42299

Installed Product: 58696268

UDI:

	Min	Max	UoM	Measurement	Result	Comments
PM section 4. Safety						
4. Safety						
Prevent electro static discharge					PASS	
PM section 6. Workflow						
On-site preparations						
Interview the representative of the site					YES	
PM section 7. System level						
System level						
Do PIQT (Remote / On-site)					PASS	
Check system parts visually					PASS	
Clean the computers Internally					PASS	
Check the dust filter					PASS	
Clean the dust filter					PASS	
Check the dust filter of the patient ventilation					PASS	
PM - Disable alerts [≥ R10]					Not Applicable (< R10)	
Check the RF enclosure					PASS	
Check the RF coils for damage					PASS	
Relative humidity Examination room	40	70	%	66.7	Pass	
Relative humidity Technical room	30	70	%	68.4	Pass	
Temperature Examination room (> R5)	18	22	°C	22	Pass	

Min	Max	UoM	Measurement	Result	Comments
-----	-----	-----	-------------	--------	----------

PM section 7. System level

System level

Temperature Technical room	15	24	°C	19.9	Pass
Check the respiratory sensor and hose					PASS
Check the version of the SPT spec file (Remote / On-site)					PASS
Compare and download the SPT spec file					PASS

PM section 8. Mains distribution Ingenia

Earth bonding examination room

SFB - Magnet	50	mOhm	11	Pass
SFB or Magnet - Patient table	100	mOhm	6	Pass
SFB or Magnet: MEU	150	mOhm	10	Pass
SFB or Magnet: PFEI	100	mOhm	4	Pass

Earth bonding tech and control room

MDU - SFB:	30	mOhm	12	Pass
MDU - RF enclosure:	30	mOhm	12	Pass
MDU - MDU (door):	50	mOhm	10	Pass
MDU or SFB: LCC cover electronic box	50	mOhm	38	Pass
MDU or SFB - Cryo compressor in LCC	100	mOhm	36	Pass
MDU or SFB - LCC cabinet	100	mOhm	40	Pass
MDU or SFB - Mains inlet unit	100	mOhm	17	Pass
MDU or SFB - DACC cabinet	100	mOhm	28	Pass
MDU or SFB - CDAS	100	mOhm	20	Pass
MDU or SFB - Reconstructor chassis	150	mOhm	24	Pass
MDU or SFB - Front cover DACC	150	mOhm	28	Pass
MDU or SFB - HOS	100	mOhm	19	Pass
MDU or SFB - Mains inlet unit ACC 1H RF amp(s)	50	mOhm	15	Pass
MDU - Mains Inlet Unit DACC	100	mOhm	17	Pass
MDU or SFB - RF amp. 2	100	mOhm	20	Pass
MDU or SFB - Circulator	100	mOhm	14	Pass
MDU or SFB - Gradient amplifier cabinet	100	mOhm	13	Pass
MDU or SFB - Front cover grad amp cabinet	150	mOhm	13	Pass
SFB or magnet - PE of SACU	100	mOhm	30	Pass
SFB or Magnet: Host/Host-recon PC chassis	100	mOhm	13	Pass

Earth Leakage Current

Leakage current (system input)	20	mA	2.12	Pass
Leakage current (transf. output)	20	mA	2.59	Pass

Min	Max	UoM	Measurement	Result	Comments
-----	-----	-----	-------------	--------	----------

PM section 9. Gradient System

Gradient System

Check the circuit to detect overheated gradient coils				PASS	
---	--	--	--	------	--

PM section 10. Liquid Cooling

Liquid cooling

Check the filling pressure	1.4	2.4	bar	2	Pass
----------------------------	-----	-----	-----	---	------

Check the hoses with liquid cooling for leaks					PASS
---	--	--	--	--	------

Liquid cooling [LCC2] [LCC2A]

Fill the circuit of the gradient coil					PASS
---------------------------------------	--	--	--	--	------

Fill the circuit of the gradient amplifier					PASS
--	--	--	--	--	------

Check the flow in the gradient coil circuit	18	22	l/minute	22	Pass
---	----	----	----------	----	------

Check the flow in the gradient amplifier	30	40	l/minute	35	Pass
--	----	----	----------	----	------

PM section 11. Refrigerator system

Refrigerator system

Check the compressor pressure F50 (dynamic)	1.9	2.2	MPa	2.2	Pass
---	-----	-----	-----	-----	------

Check the cryogenic performance					PASS
---------------------------------	--	--	--	--	------

PM section 12. Magnet system

Magnet system

Check the labels and covers of the ERDU buttons					PASS
---	--	--	--	--	------

Check the outlet of the helium venting system					PASS
---	--	--	--	--	------

Check the service history of the coldhead (48 months)					≤ 48 Months Ago
---	--	--	--	--	-----------------

Record the last date coldhead service was done					Nov 11, 2022
--	--	--	--	--	--------------

Check the ERDU with the magnet ON					PASS
-----------------------------------	--	--	--	--	------

Check Differential Pressure Switch (DPS)					PASS
--	--	--	--	--	------

PM section 13. Patient support

Patient support [IMT]

Check the patient alarm					PASS
-------------------------	--	--	--	--	------

Check the patient comfort lights					PASS
----------------------------------	--	--	--	--	------

Clean the DCI RF coil sockets and connector					PASS
---	--	--	--	--	------

Clean the PSU					PASS
---------------	--	--	--	--	------

Min	Max	UoM	Measurement	Result	Comments
-----	-----	-----	-------------	--------	----------

PM section 14. System Meas. & Adjustments

System Meas. & Adjustments

Check the body coil transmit characteristics				Not applicable	
--	--	--	--	----------------	--

Run the software tests for PM

System selftest				PASS	
PFEI peripheral test				PASS	
Patup self test				PASS	
PSU2 selftest (PSU(2) fan speed test)				PASS	
TX att. calibration				PASS	
Max kW calibration body coil				PASS	
RF amp TX1 PMU test body coil				PASS	
RF amp TX2 PMU test body coil				PASS	
RF amp controlloop calibration				PASS	
F0 determination				PASS	
RF power ref. cal. body coil				PASS	
Pickup coil tripl.cal.body coil				PASS	
Multitransmit RF calibration				PASS	
Max kW calibration TR-Knee coil				PASS	
RF power ref. cal. TR-Knee coil				PASS	
Pickup coil tr.cal TR-Knee coil				PASS	

System measurements and adjustments

Check the VSWR (3.0T multi-transmit TX1)	5		3.16	Pass	
Check the VSWR (3.0T multi-transmit TX2)	5		2.84	Pass	
Use the analyzer				PASS	
Run the RF spurious test				PASS	
Run the spurious noise test batch [$\leq$ R5.5]				Not Applicable $\geq$ R5.6	
Analyze the result of the spurious noise test batch [$\leq$ R5.5]				Not Applicable $\geq$ R5.6	

PIQT

Do PIQT				PASS	
Generate a report with PIQT results				PASS	

PM section 15. Operator Console & Viewing

Operator Console & Viewing (All)

Clean the mouse, keyboard, monitor of host and EWS				PASS	
Check the audio communications				PASS	

Min	Max	UoM	Measurement	Result	Comments
-----	-----	-----	-------------	--------	----------

PM section 16. Software

Software

Check and update the EAS software (e-Alert) from Incenter to the latest version (select PASS or FAIL if micro SD Card installed)	e-ALERT NOT INSTALLED
Clean up disks	PASS
Defragment the hard disks [Windows XP]	Not Applicable (Not Windows XP)
Make a system back-up	PASS
Export the configuration	PASS

PM section 19. PM end

19. PM end

Handover	PASS
Administration	PASS

Philips Representative Signature:

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Engineer Name: Takuya Nishio
Date: Aug 2, 2024

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